

CONTROL CENTER CARE BED THE PERFECT COMBINATION

#16444 V1.1 EN



BED CONTROL



RemoBED

The **RemoBED Solutions** are suitable for many care beds, especially for the care bed „Dali“ from Burmeier.



csslabs.de/remobed

&

SECURE PAGING



AssistX

Assured calling with the **AssistX Call System** in intensive care, at home, in residential groups, in nursing homes and clinics.



csslabs.de/assistx



BED CONTROL

FOR ALL!



CSS MicroSystems GmbH
In den Hofäckern 16
74189 Weinsberg
Germany
www.csslabs.de



RemoBED Assistive Bed Control

Adjusting the care bed independently:
By pressing buttons, by speaking
or with eye & head control

RemoBED

Assistive Bed Control

BED CONTROL FOR ALL!

Being able to adjust the care bed independently is a big step to more quality of life.

RemoBED solutions help everyone who cannot operate the standard handset.

RemoBED solutions are simple and flexible to mount to the care bed.

RemoBED solutions at a glance:



RemoBED Functions

With **RemoBIF** these functions are adjustable (up and down):

The backrest

The footrest

The bed height

(Picture: The „Dali“ care bed from Burmeier, but also available for many other care beds)



RemoBOW

By pressing buttons

With **RemoBOW** it is possible to operate the bed with big buttons. The backrest and footrest can be adjusted by simply pressing these buttons. An **AssistX Call System** and switching a lamp on and off can be integrated.



RemoSPEAK

By speech control

With **RemoSPEAK** it is possible to control the bed with simple voice commands. The backrest and footrest as well as the bed height can be adjusted. An **AssistX Call System** and switching a lamp on and off can be integrated.



RemoBIF

With eye & head control

Using a communication device, it is possible to operate the care bed with eye or head movements. The backrest and footrest as well as the bed height can be adjusted.



RemoBIF und BIF8+
Universal bed control interface
for environment control systems